



Silicone Foam Border Dressing & Silicone Foam Border Sacral Dressing



- Highly absorbent primary dressing
- Maintains balanced moist wound environment
- Gentle silicone wound face & border
- Soft & conformable
- Can be used on a variety of wound types
- Valuable tool in wound management

PIVOT Silicone Foam Border Dressing is a gentle, conformable foam dressing with a wound facing perforated silicone layer and silicone adhesive border. The customized foam matrix maximizes exudate absorption and retention while maintaining structure and conformability. The perforated silicone layer is gentle and does not adhere to the wound thus facilitating easy repositioning or removal. **PIVOT** Silicone Foam Border Dressing helps create and maintain a moist wound environment optimal for wound healing.

ITEM DESCRIPTION	ITEM NUMBER	SIZE: INCH (cm)	UNITS/CARTON	HCPCS
Silicone Foam Border Dressing	SF2X2B	2"x2" (5cmx5cm) Pad Size: 1"x1" (2.5cmx2.5cm)	10 Dressings	A6212
Silicone Foam Border Dressing	SF3X3B	3"x3" (7.5cmx7.5cm) Pad Size: 1.5"x1.5" (3.8cmx3.8cm)	10 Dressings	A6212
Silicone Foam Border Dressing	SF4X4B	4"x4" (10cmx10cm) Pad Size: 2.5"x2.5" (6.2cmx6.2cm)	10 Dressings	A6212
Silicone Foam Border Dressing	SF6X6B	6"x6" (15cmx15cm) Pad Size: 4.2"x4.2" (10.5cmx10.5cm)	10 Dressings	A6213
Silicone Foam Border Sacral Dressing	SF7X7B	7"x6.8" (17.75cmx17cm) Pad Size: 5"x4.8" (12.5cmx12cm)	10 Dressings	A6213