



## Silver Alginate Dressing



- **Metallic silver particles inhibit the growth of Gram positive and Gram negative bacteria within the dressing material**
- **Highly absorbent up to 20 times its weight**
- **All-natural**
- **Maintains a moist wound environment**
- **Adaptable to wound size & shape**
- **Easy to apply, reposition, and remove**
- **Valuable tool in wound management**

**PIVOT** Silver Alginate Dressing is a fibrous, highly absorbent, all-natural biodegradable material derived from seaweed. It is highly absorbent and absorbs up to 20 times its weight. The dressing absorbs exudates and forms a gel-like covering over the wound, which maintains a moist environment for optimal wound healing. The gel-like surface also keeps the dressing from adhering to the wound, which makes dressing removal painless.

Metallic silver particles are added to the Calcium Alginate fibers that inhibit the growth of bacteria in the dressing. Each dressing contains 50ug/100cm<sup>2</sup> metallic silver particles.

ITEM DESCRIPTION	ITEM NUMBER	SIZE: INCH (cm)	UNITS/CARTON	HCPCS
Silver Alginate Dressing	SA2X2	2"x2" (5cmx5cm)	10 Dressings	A6196
Silver Alginate Dressing	SA4X4	4"x4" (10cmx10cm)	10 Dressings	A6196
Silver Alginate Dressing	SA4X5	4"x5" (10cmx12.5cm)	10 Dressings	A6197
Silver Alginate Dressing	SA6X6	6"x6" (15cmx15cm)	10 Dressings	A6197
Silver Alginate Rope Dressing	SA1X12	1"x12" (2.5cmx30cm)	10 Dressings	A6199